

Test Kit Instructions

January 03, 2020

VICAM AFLATEST

QUANTITATIVE AFLATOXIN TEST

FORWARD

The instructions presented in this document cover only the procedure for performing the analytical test for official inspections. For questions regarding this procedure, contact Dr. Ajit Ghosh of the Technology and Science Division (TSD) by phone at 816-891-0417 or email at Ajit.K.Ghosh@usda.gov.

Refer to the Mycotoxin Handbook for information on use of this test kit in official FGIS inspections including sampling, general sample preparation, reporting and certification of test results, laboratory safety, and hazardous waste management. For questions regarding these policies and/or instructions, contact Patrick McCluskey of the Policies, Procedures, and Market Analysis Branch (PPMAB) by phone at 816-659-8403 or email at Patrick.J.McCluskey@usda.gov.

The U.S. Department of Agriculture (USDA) prohibits discrimination in its programs on the basis of race, color, national origin, sex, religion, age, disability, political beliefs, and marital or familial status. (Not all prohibited bases apply to all programs.) Persons with disabilities who require alternate means for communication of program information (Braille, large print, audiotape, etc.) should contact USDA's TARGET Center at (202) 720-2600 (voice and TDD).

To file a complaint, write to the USDA, Office of Civil Rights, Room 326-W, 1400 Independence Avenue, SW, Washington, DC 20250-9410, or call (202) 720-5964 (voice and TDD). USDA is an equal employment opportunity employer.

Contents

1.	GENERAL INFORMATION	3
2.	PREPARATION OF TESTING MATERIALS AND EQUIPMENT	4
3.	EXTRACTION PROCEDURES	7
4.	TEST PROCEDURES	9
5.	REPORTING AND CERTIFYING TEST RESULTS	11
6.	STORAGE CONDITIONS AND PRECAUTIONS	11
7.	EQUIPMENT AND SUPPLIES	12
8.	REVISION HISTORY	14

1. GENERAL INFORMATION

The VICAM AflaTest method uses immunoaffinity column chromatography to isolate aflatoxins B1, B2, G1, and G2 for fluorometric detection and quantification at the parts-per-billion (ppb) level.

Approved Test Kit Information	
Test Kit Vendor:	VICAM, A Waters Business (508) 482-4935
Test Kit Name:	VICAM AflaTest
Product Number:	12022
Effective Date of Instructions:	12/18/2019
Type of Service:	Quantitative
Conformance Range:	5.0 – 300 ppb
Number of Analyses to Cover Conformance Range:	1
Approved Commodities:	Corn (including dent or field corn, corn meal, corn flour, cracked corn, corn grits or polenta, and corn screenings), corn gluten meal, distillers dried grain with solubles (DDGS), milled rice (including brewer's rice and glutinous rice), popcorn, and sorghum
Extraction method:	Blend 50 grams (g) sample with 5 g of sodium chloride and 100 milliliters (mL) methanol:water (80:20, v:v) for 1 minute (for DDGS use 10 g sodium chloride and 200 mL methanol:water (80:20, v:v). For corn gluten meal, blend 50 g sample with 10 g sodium chloride and 200 mL methanol:water (60:40, v:v) for 3 minutes.
Test Format:	Immunoaffinity column with fluorometric detection
Detection Method:	VICAM Series 4 EX fluorometer

2. PREPARATION OF TESTING MATERIALS AND EQUIPMENT

Before beginning the procedures, please read the “Precautions Section” on page 11.

a. VICAM 4EX Fluorometer Calibration.

(1) General.

- (a) To ensure accurate results, calibrate the fluorometer before use each day and verify calibration using the yellow vial.
- (b) Turn the fluorometer on with the On/Off switch located on the rear panel. Once the fluorometer is turned on, it may be left on until close of business for the day.
- (c) After turning the fluorometer on, it will identify itself and perform a set of self-tests. If any error message appears, consult the operator’s manual or call Vicam at 1-800-338-4381.

b. Calibration Procedures.

- (1) Set the date and time by pressing the OPTIONS key until “Settings” appears on the display; then press ENTER. “Set Date and Time” will appear on the display. Press ENTER, and follow the prompts.
- (2) Press the STOP key to return the fluorometer to the “VICAM ready” mode.
- (3) Press the OPTIONS key until “Calibrate Test” appears on the display. Press ENTER. Press the SELECT TEST key until “AflaTest FGIS” is on the display. Press ENTER, and follow the prompts.
- (4) When prompted to insert a calibration vial, clean the vial with a Kim wipe and insert it into the bottom of the well. Ensure the vial is fully inserted and touches the bottom of the well.
- (5) Enter the correct calibration value (see table below) for the high calibrator (red vial) and low calibrator (green vial).

Calibration Vials	Series 4 EX
Red	160
Green	-1.0
Yellow	79 ± 5

- (6) Check the calibration by testing the yellow vial.
- (7) Press the SELECT TEST key. The display should read "AflaTest FGIS." Press ENTER, and then insert the yellow vial. Record the results for the yellow vial on the work record.
- (8) If the value of the yellow calibration vial is not within FGIS specifications, repeat the calibration process (steps b through f); then check the yellow vial again. If the reading for the yellow vial reads less than 74 or above 84, contact Vicam and let TSD know about the incident.
- (9) Calibrate the fluorometer at least once a day or before analyzing any test samples if more than 1 day has elapsed since the last test using the yellow vial.

c. Reagent Testing.

After calibrating the fluorometer, check the following reagents for background fluorescence by performing the procedures listed below.

- (1) Methanol: Dispense 2.0 milliliters (mL) of HPLC grade methanol into a clean cuvette and place the cuvette in the fluorometer. The digital display should read zero (0).
- (2) Water: Dispense 2.0 mL of deionized or distilled water into a clean cuvette and place the cuvette in the fluorometer. The digital display should read zero (0).
- (3) If any of the readings is greater than zero (0), replace reagent from a pure source, recheck fluorometer calibration, or replace cuvette.

d. Developer Solution.

- (1) Using a 5 mL pipettor, add 5 mL of AflaTest concentrated developer solution to 45 mL of distilled or deionized water (measured using a 50 mL graduated cylinder). Mix well. The concentrated developer solution should have a slight yellowish-

brown color. If the solution is colorless, it has lost potency and should not be used. When preparing smaller quantities of the solution be sure to maintain the one part concentrated developer to nine parts distilled or deionized water ratio.

- (2) Measure 1.0 mL of HPLC grade methanol and 1.0 mL of dilute developer solution with a 1 mL pipettor and combine in a clean VICAM borosilicate glass cuvette. Place the cuvette in the fluorometer to check for background fluorescence. The digital display reading should read zero (0). If the reading is greater than 0, check each reagent separately to determine which reagent is causing the problem and replace it, recheck calibration, or replace cuvette.
- (3) Label the dilute developer solution bottle with the date and time of preparation and the technician's initials. Fresh developer solution must be prepared daily and used within 8 hours.
- (4) Label each stock bottle of concentrated developer with the date it was first opened. Do not use concentrated developer after 30 days have elapsed.

e. Extraction Solvent: Methanol:Water (80:20,v:v).

Premixed methanol:water (80:20:v:v) solution may be purchased from VICAM (part # 100000211-20L or part # 100000212-4L).

One may prepare extraction solvent in a well-ventilated area as follows:

- (a) Using a 1000 mL graduated cylinder, measure 800 mL of methanol (ACS grade or better), and pour it into a clean carboy with spigot.
- (b) Using a 250 mL graduated cylinder, measure 200 mL deionized or distilled water, and add to the methanol. Shake until completely mixed. Allow to reach room temperature before use.
- (c) Label the container stating the mixture (methanol:water (80:20 v:v)), date of preparation, and technician's initials.
- (d) Store this solution at room temperature in a tightly closed container until needed.

f. Extraction Solvent: Methanol:Water (60:40 v:v) for corn gluten meal only.

Premixed methanol:water (60:40, v:v) may be purchased from VICAM (part # 100000213-20L or part # 100000214-4L).

One may prepare it as follows.

- (a) Using a 1000 mL graduated cylinder, measure 600 mL of methanol (ACS grade or better), and pour it into a clean carboy with spigot.
- (b) Using a 500 mL graduated cylinder, measure 400 mL deionized or distilled water and add to the methanol. Shake until completely mixed.
- (c) Label the container stating the mixture (methanol:water (60:40, v:v)), date of preparation, and the technician's initials.
- (d) Store this solution at room temperature in a tightly closed container until needed.

g. Column chromatography setup.

- (1) The Test Procedures will require a VICAM pump stand with a glass syringe barrel reservoir and a manually operated syringe pump or an adjustable automated air pump to push the sample extract and reagents through the AflaTest column.
- (2) Before beginning Test Procedures, remove the top and bottom caps from the AflaTest affinity column and gently shake the buffer solution from the top of the column.
- (3) When attaching the column to the glass syringe barrel on the pump stand, use the appropriate column coupling. Create a coupling by snipping off the bottom of the top cap after removing it from the column. The coupling can be reused.

3. **EXTRACTION PROCEDURES**

Place 50 grams (± 0.2 grams) of ground sample into blender jar and add the appropriate reagents before blending, filtering and diluting (step b).

- a. Reagents for corn (including dent or field corn, corn meal, corn flour, cracked corn, corn grits or polenta, and corn screenings), milled rice (including brewer's rice and glutinous rice), popcorn, and sorghum:

- (1) Add 5 grams (g) of analytical, USP grade sodium chloride (NaCl) or food grade non-iodized salt.
 - (2) Using a 250 mL graduated cylinder (VICAM # 20250), add 100 mL of the methanol:water (80:20, v:v) extraction solvent.
- b. Reagents for distillers dried grains with solubles:
- (1) Add 10 g of analytical, USP grade sodium chloride (NaCl) or food grade non-iodized salt.
 - (2) Using a 250 mL graduated cylinder, add 200 mL of the methanol:water (80:20, v:v) extraction solvent.
- c. Reagents for corn gluten meal:
- (1) Add 10 g of analytical, USP grade sodium chloride (NaCl) or food grade non-iodized salt.
 - (2) Using a 250 mL graduated cylinder, add 200 mL of the methanol:water (60:40, v:v) extraction solvent.
- d. Blending, Filtering and Dilution
- (1) Cover the jar and blend at high speed for 1 minute. (for corn gluten meal blend for 3 minutes)
 - (2) Pour the extract through a clean funnel lined with fluted filter paper into a clean VICAM disposable beaker labeled with sample identification.
 - (3) After collecting approximately 20 mL of filtered extract, carefully dispose of filter and its contents. Filtered extract must be used within 30 minutes.
 - (4) Pipette 10 mL of distilled or deionized water using 5 mL pipettor (VICAM # G4055) into an extraction beaker (VICAM # 600000827).
 - (5) Using a new pipette tip, add 5 mL of filtered extract and mix it by vortex.
 - (6) Line a small, clean funnel with a VICAM microfiber filter. Fold filter gently without making a sharp crease to avoid breaking it. Pour diluted extract through the lined funnel into a 12 x 75 mm disposable borosilicate glass cuvette (VICAM # 34000) or hard plastic beaker (VICAM # 36010). This is the diluted filtered extract and ready for the analysis.

- (7) Immediately carry out the test procedure.

4. TEST PROCEDURES

Perform the appropriate sample cleanup procedures before proceeding to step d (Separation and Measurement) on page 10.

- a. Sample cleanup for corn (including dent or field corn, corn meal, corn flour, cracked corn, corn grits or polenta, and corn screenings), milled rice (including brewer's rice and glutinous rice), popcorn and sorghum:

- (1) Using 1 mL pipettor, pipette 1.0 mL of the diluted filtered extract into the headspace of the AflaTest column.
- (2) Use the column coupling to attach the AflaTest column to the bottom of the glass syringe barrel on the pump stand. Place a waste collection vessel under the column outlet. Insert the coupling on the end of pump tubing into the glass syringe barrel and activate the pump. Maintain a steady flow rate of approximately 1 drop of extract per second until all the extract has passed through the column. Maintain pressure on the column until air passes through the column.
- (3) Detach the test column from the syringe barrel and fill the column headspace with deionized or distilled water (about 1 mL). Reattach column to the syringe barrel and activate the pump. Maintain a steady flow rate of 1 to 2 drops of water per second until all the water has passed through the column.
- (4) Repeat step 3. Maintain pressure on the column until air passes through the column.

- b. Sample cleanup for distillers dried grains with solubles:

- (1) Pipette 1.0 mL of the diluted filtered extract into the headspace of the AflaTest column.
- (2) Use the column coupling to attach column to the bottom of the glass syringe barrel on the pump stand. Place a waste collection vessel under the column outlet. Insert the coupling on the end of the pump tubing into the top of the glass syringe barrel, and activate the pump. Maintain a steady flow rate of approximately 1

drop of extract per second until all the extract has passed through the column.

- (3) Detach the column from the glass syringe barrel. Repeat steps 1 and 2. Maintain pressure on the column until air passes through the column
- (4) Detach the column from the syringe barrel and fill the column headspace with deionized or distilled water (about 1mL). Reattach the column to the syringe barrel, add 5 mL water to the syringe barrel and activate the pump. Maintain a steady flow rate of 1 to 2 drops of water per second until all the water (6 ml total) has passed through the column. Maintain pressure on the column until air passes through the column.

c. Sample Cleanup for Corn Gluten Meal:

- (1) Use the column coupling to attach the column to the glass syringe barrel on the pump stand. Place a waste collection vessel under the column outlet. Add 4 mL of the diluted filtered extract to the syringe barrel. Insert the coupling at the end of the pump tubing into the top of the glass syringe barrel and activate the pump. Maintain a steady flow rate of approximately 1 drop of extract per second until all the extract has passed through the column.
- (2) Add 10 mL of deionized or distilled water to the syringe barrel and activate the pump. Maintain a steady flow rate of 1 to 2 drops of water per second until all the water has passed through the column.
- (3) Repeat step 2.

d. Separation and Measurement.

- (1) Place a clean VICAM borosilicate glass cuvette under the outlet of the IAC column.
- (2) Detach the column and using 1 mL pipettor, pipette 1 mL of HPLC grade methanol directly into the column headspace. Reattach the column and activate the pump. Apply just enough air pressure to pass the methanol through the column at a rate of approximately 1 drop **every 2 seconds**. Collect all the methanol (1mL) in the cuvette.

- (3) Using 1 mL pipettor, add 1.0 mL of freshly made developer solution directly to the sample eluate in the cuvette. Vortex for 3 to 5 seconds. Clean the outside of the cuvette quickly with a Kim wipe.
- (4) Immediately place the cuvette in a previously calibrated fluorometer. Reader reads sample concentration after 60 seconds.
- (5) Sample results above 300 ppb are reported as “Aflatoxin exceeds 300 ppb.”

5. REPORTING AND CERTIFYING TEST RESULTS

Refer to the Mycotoxin Handbook for reporting and certification of test results. For questions regarding these instructions, contact Patrick McCluskey (816-659-8403 or Patrick.J.McCluskey@usda.gov).

6. STORAGE CONDITIONS AND PRECAUTIONS

a. Storage conditions.

- (1) Store AflaTest columns at room temperature (64° to 77° F).
- (2) The standard solutions in the three standard vials (red, green, and yellow) degrade slowly in the presence of light. Keep vials in closed case. About 2 months before the expiration date, obtain a new set of AflaTest FGIS Calibration Standards (VICAM product # 33030).
- (3) Store the developer concentrate in a tightly closed bottle in a cool, dry area and away from sunlight.

b. Precautions.

- (1) The sample extract passed over the AflaTest column must be clear. If the diluted filtrate turns cloudy, filter again using a new glass microfiber filter before proceeding with testing.
- (2) Use the recommended equipment for procedures:
 - (a) Collection of extracts and filtrates: VICAM disposable hard plastic beakers (VICAM # 36010), Extraction tubes (VICAM # 600000827) or clean glass

containers. Soft plastic may leach fluorescent materials.

- (b) Fluorometer cuvettes: VICAM 12 x 75 mm borosilicate glass cuvettes (VICAM # 34000). Note: Cuvettes from other sources may emit significant fluorescence, resulting in inaccurate readings. Keep cuvette's optical surface free of lint, fingerprints, etc.
- (3) Add diluted developer solution directly to the cuvette containing the eluate from the AflaTest column. Use freshly made (within 8 hours) solution only. Mix well by vortexing or flicking bottom of the tube with fingers. DO NOT shake or invert cuvette using a thumb or finger as a cap. Doing this may elute fluorescent materials.
- (4) Carry out extraction of sample in a well-ventilated area.

7. EQUIPMENT AND SUPPLIES

a. Materials required.

- (1) VICAM Series 4EX Fluorometer (VICAM # FLSEREX)
- (2) AflaTest FGIS Calibration Standards (VICAM # 33030)
- (3) Pump stand with air pump (VICAM # 21020, 21030, 21040, 21045, G4061 or G1104)
- (4) AflaTest developer solution (VICAM # 32010)
- (5) AflaTest immunoaffinity columns (VICAM # 12022)
- (6) Distilled or deionized water
- (7) ACS grade or better methanol
- (8) HPLC grade methanol, 4L (VICAM # 100000219)
- (9) Premixed methanol:water (80:20, v:v) solution (VICAM # 100000211- 20L or #100000212-4L)
- (10) Premixed methanol:water (60:40, v:v) solution for corn gluten meal only (VICAM # 100000213- 20L or #100000214-4L)

- (11) USP grade sodium chloride (NaCl) or food grade non-iodized salt (VICAM # G1124)
- (12) Automatic pipettor
 - (a) 1 mL capacity for methanol (VICAM # 20501)
 - (b) 1 mL capacity for developer (VICAM # 20600)
- (13) Blender base (VICAM # 20201)
- (14) Eberbach Glass Blender Jar, 500mL (VICAM # 20300)
- (15) Cuvettes, disposable 12 x 75 mm borosilicate glass tube (VICAM # 34000)
- (16) Cuvette rack (VICAM # 21010)
- (17) Disposable beakers (VICAM # 36010)
- (18) Extraction tubes (40mL) pack of 25 (VICAM # 600000827)
- (19) Fluted filter paper (VICAM # 31240)
- (20) Funnel for fluted filter paper (VICAM # 36022)
- (21) Funnel for glass microfiber filter paper (VICAM # 36020)
- (22) Glass microfiber filter paper (VICAM # 31955)
- (23) Graduated cylinder, 50 mL glass (VICAM # 20050)
- (24) Graduated cylinder, 250 mL plastic (VICAM # 20250)
- (25) Graduated cylinder, 500 mL plastic (VWR# 470188-960)
- (26) Graduated cylinder, 1000 mL plastic (VWR # 65000-012)
- (27) Syringe, glass 10 mL (VICAM # 34010)
- (28) Water wash bottle (VICAM # 20700)
- (29) Micro-pipettor, 1mL (VICAM # G4033)
- (30) Micro-pipette Tips for 1mL Micro-pipettor, 100 (VICAM # 20656)
- (31) Adjustable pipettor, 5mL (VICAM # G4055)

- (32) Finntips for 5mL adjustable pipettor (VICAM # G4060)
- (33) Kim wipes (VICAM # 31967)
- (34) A package of all the basic equipment can be purchased from VICAM (VICAM # G8001)

8. REVISION HISTORY

Effective: 01/03/2020

- Corn gluten meal, distillers dried grain with solubles (DDGS), milled rice (including brewer's rice and glutinous rice), popcorn, and sorghum were approved as the additional commodities. The test procedure has been updated to accommodate these additional commodities.

Effective: 11/13/2019